

KingScope™ 100 Brinell Microscope

METAL HARDNESS TESTING MADE EASY

The KingScope™ 100 is the most accurate manual microscope on the market with resolution at .01 mm. With its ergonomic design and small nose piece, the KingScope™ 100 can be used in small areas. Its lightweight design is extremely rugged with a military drop-tested case to accommodate harsh conditions.

Focusing the KingScope™ 100 is fast and easy, and the built-in LED provides additional light when needed



ACCURATE

- 20x total magnification
- Accurate to nearest 0.01 mm
- High-quality optics
- Meets ASTM E-10 Type A Scope requirements



VERSATILE

- Narrow nose piece fits easily into small areas
- Lightweight yet rugged ergonomic design for easy handling
- LED for easy reading in all light conditions
- Suitable for both laboratory and shop
- Military drop-tested case and penlight to accommodate for harsh conditions



Microscope stand for convenient and easy access. The stand can attach to any of our portable Brinell tester bases or be used as a stand alone.