

## Brinell KingScope™ Scope

## METAL HARDNESS TESTING MADE EASY

Constructed from stainless steel, the rugged and optically reliable KingScope is the most versatile scope on the market today. Featuring a 20x pre-focused lens, the KingScope has a narrow nose piece which easily fits into tight recesses, resulting in less grinding on castings, billets and dies. For added stability when performing flat work, a slip-on base adapter is included. A side opening in the scope allows plenty of natural light for viewing, and a cordless movable pen light can be used in dim conditions.

The KingScope meets ASTM E-10 specifications and the calibration is traceable to NIST standards.



## **ACCURATE**

- 20x total magnification
- Accurate to nearest 0.1mm
- Pre-focused lens allows for a fast and accurate reading
- High-quality optics
- Meets ASTM E-10 standards



## **VERSATILE**

- Narrow nose piece fits easily into small areas
- Portable design
- Side opening and cordless light allows for easy reading in all light conditions
- Suitable for both laboratory and shop
- Sturdy stainless steel construction



